




Determination of pK_a value of a weak acid


Aim of the experiment:

To determine the pK_a value of a weak acid potentiometrically using a pH meter



https://www.bing.com/images/search?view=detailV2&ccid=67X5U41t&id=782B4A43F25B5A14E9D075393293752E80CB0D59&thid=OIP.67X5U41tzFXUESZ5-nEWowHaHa&mediaurl=https%3a%2f%2fimg.medicaexpo.com%2fimages_me%2fphoto-g%2f100041-16315454.jpg&exph=1500&expw=1500&q=Ph+Metre&simid=608013421153759747&FORM=IRPRST&ck=C48F69EE00DB825B1440EB26C86DE87F&selectedIndex=0&itb=0

pH indicator electrode - Combined glass electrode



Scheme of typical pH glass electrode

1. a sensing part of electrode,
2. a bulb made from a specific glass
sometimes electrode contain small amount
of AgCl precipitate inside the glass
electrode
- 3 internal solution, usually 0.1M HCl for pH
electrodes
- 4.internal electrode, usually silver chloride
electrode or calomel electrode

- 5.body of electrode, made from non-
conductive glass or plastics.
- 6.reference electrode, usually the same type
as 4
- 7.junction with studied solution, usually
made from ceramics or capillary with
asbestos or quartz fiber.

Principle and outline:

A weak acid, like acetic acid dissociates in its aqueous solution only slightly, as follows;



The equilibrium constant of dissociation of the acid known as dissociation constant of the acid, K_a is given by:

$$K_a = \frac{[\text{H}_3\text{O}^+][\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]}$$

The solution is considered to be very dilute and activity of water in very dilute solutions is assumed to be unity in the above equation.

$$\therefore [\text{H}_3\text{O}^+] = \frac{K_a [\text{CH}_3\text{COOH}]}{[\text{CH}_3\text{COO}^-]}$$

From the above equation, taking $-\log$ on both sides:

$$-\log [\text{H}_3\text{O}^+] = -\log K_a - \log \frac{[\text{CH}_3\text{COOH}]}{[\text{CH}_3\text{COO}^-]}$$

$$\text{i.e. pH} = \text{p}K_a + \log \frac{[\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]}$$

This is Henderson's equation. When acetic acid (titrate) is titrated with NaOH solution (titrant), pH of the titrated solution increases in accordance with the above equation; since concentration of sodium acetate, $[\text{CH}_3\text{COO}^-]$ increases and concentration of unreacted acetic acid, $[\text{CH}_3\text{COOH}]$ decreases as the titration progresses. When the acid is half - titrated,

$$[\text{CH}_3\text{COO}^-] = [\text{CH}_3\text{COOH}]$$

and hence $\text{pH} \equiv \text{p}K_a$

That is, $\text{p}K_a$ of the acid is pH of the half-titrated acid solution or pH of the titrated solution at half way to the equivalence point or end-point.

Materials Needed:

- pH meter
- Combined pH electrode
- Buffer solutions of known pH (e.g., pH 4)
- Beakers
- Weak acid solution (50 mL)
- Sodium hydroxide (NaOH) solution
- Semi-micro burette

Procedure:



I. Standardization of pH Meter - Dipping in Buffer:

- Dip the combined electrode into a beaker containing the buffer solution (e.g., pH 4).
- Connect the electrode to the pH meter.
- Observe the reading on the pH meter. It should correspond to the pH of the buffer solution (pH 4).
- If the reading does not match the known pH of the buffer, adjust the control knob on the pH meter until the display shows the correct pH value.
- Once standardized, the pH meter is ready for use.

II. Potentiometric titration

- Pipette 50 mL of the weak acid into a clean beaker.
- Rinse the combined electrode with distilled water and then dip it into the weak acid solution.
- Connect the electrode to the pH meter.
- Record the initial pH of the weak acid solution.
- Using a semi-micro burette, add exactly 1 mL of NaOH to the acid solution.
- Stir the solution gently to ensure proper mixing.
- Measure and record the pH after each addition of NaOH.
- Continue adding 1 mL increments of NaOH and recording the pH until a significant change in pH is observed.
- Continue the titration up to a total of 10 mL of NaOH.
- After reaching 10 mL, freshly pipette out another 50 mL of the weak acid.
- Record the pH after the addition of smaller increments (0.1 mL) of NaOH in the region where large changes in pH were previously observed.
- Continue this process until the change in pH becomes small again.
- Plot the data on a graph with the volume of NaOH on the X-axis and pH on the Y-axis.
- Analyze the graph to identify the inflection point, which corresponds to the pK_a of the weak acid.

Graph:



pKa value of acetic acid =

Applications of pH meter

1. Agriculture

- Soil testing for acidity/alkalinity.
- Optimizing fertilizer use.

2. Environmental Monitoring

- Assessing water quality in ecosystems.
- Detecting pollution levels.

3. Food and Beverage Industry

- Ensuring quality and safety in food products.
- Monitoring fermentation processes in brewing and winemaking.

4. Healthcare and Pharmaceuticals

- Measuring blood and urine pH for diagnostics.
- Controlling pH in drug formulation.

Instrument Calibration and Measurement of pH in water
by pen type pH Meter



Thank you

